



PeliKine Compactä Human IFN γ ELISA kit

288 tests

An enzyme immunoassay for the quantitative
determination of human interferon gamma

PRODUCT INFORMATION

Central Laboratory of the
Netherlands Red Cross
Blood transfusion Service
PO box 9190
1006 AD Amsterdam1006
AD Amsterdam

Cat.No. M1933

For the Netherlands:

Tel. 020.512 3599

Fax. 020.512 3570

For other countries:

contact your local distributor.

E-mail: clbreagents@clb.nl

Website: <http://www.CLB.nl>

Protocol summary and checklist PeliKine compactä human IFN γ ELISA kit

Day 0:

- Bring coating antibody to room temperature (18-25 °C).
- Prepare coating buffer.
- Dilute coating antibody 1:100 in coating buffer, add 100 μ l to all wells, cover the plate(s) and incubate overnight at room temperature.

Day 1:

- Bring all reagents, with the exception of streptavidin-HRP, to room temperature.
- Prepare blocking buffer
- Wash the plate(s) five times with PBS.
- Add 200 μ l blocking buffer to all wells and incubate for one hour at room temperature.
- Prepare standard and sample dilutions.
- Prepare washing buffer.
- Wash the plate(s) five times with washing buffer.
- Leaving the substrate blank wells empty, add 100 μ l of standard and sample dilutions to the appropriate wells, cover the plate(s) and incubate for one hour at room temperature.
- Dilute biotinylated IFN γ antibody 1:100 in dilution buffer.
- Wash the plate(s) five times with washing buffer.
- Leaving the substrate blank wells empty, add 100 μ l of the diluted biotinylated IFN γ antibody to all wells, cover the plate(s) and incubate for one hour at room temperature.
- Dilute the streptavidin-HRP conjugate 1:10,000 in dilution buffer.
- Wash the plate(s) five times with washing buffer.
- Leaving the substrate blank wells empty, add 100 μ l of the streptavidin-HRP conjugate to all wells, cover plate(s) and incubate for 30 minutes at room temperature.
- Just before use, prepare substrate solution.
- Wash the plate(s) five times with washing buffer.
- Add 100 μ l substrate solution to all wells, including the substrate blank wells, and incubate for 30 minutes at room temperature in the dark.
- Add 100 μ l stop solution to all wells and read the plate at 450 nm.
- Calculate the amount of IFN γ in the samples.

	1	2	3	4	5	6	7	8	9	10	11	12
A	S1	S1										
B	S2	S2										
C	S3	S3										
D	S4	S4										
E	S5	S5										
F	S6	S6										
G	S7	S7										
H	S8	S8									B	B

Plate plan proposed for the Pelikine compact™ human IFN γ ELISA kit:

Key: B: substrate blank S1-S8: IFN γ standards 0 - 500 pg/ml Empty: samples

I. INTRODUCTION

At this moment fifteen interferon α (IFN α), one interferon β (IFN β) and one interferon γ (IFN γ) have been reported. IFN γ is produced during an immune response by CD8⁺, NK, $\gamma\delta$ and TH1 T helper cells. It differs structurally and functionally from IFN α and IFN β ; binds to distinct receptors and has pronounced immuno-regulatory effects, including activation of macrophages to enhance phagocytosis and tumour killing capability, activation and growth enhancement of cytolytic T-cells and NK-cells, and induction of class II MHC antigen and Fc γ receptor on macrophages and many other cell types.

IFN γ also regulates humeral immune responses: it induces immunoglobulin secretion by activated B-cells stimulated with IL-2 and potentiates IL-4 induced proliferation of human B-cells. IFN γ has documented antiviral and antiprotozoal activities, although IFN α and IFN β seem to have more potent antiviral activities than IFN γ .

Several substances originally described for their biological activities have been identified as IFN γ ; macrophage activating factor (MAF), T-cell replacing factor (TRF), Type II interferon and immune interferon.

Bioassays for the quantification of IFN γ , based on cytopathic reductive effects of IFN γ on cultured cells have been used for several years. In this assay IFN γ reduces the killing of a target cell line such as L929 (murine), HEp2C or A549 (human) cells by for example, encephalomyocarditis virus. An alternative assay system involves measurement of induction of HLA-DR antigens on tumour cells, which can be detected in a cell ELISA. These assays, although sensitive, are time consuming and might be susceptible to interference by other substances.

The Pelikine compact™ human IFN γ ELISA kit has been developed for faster, more reproducible and specific quantification of human IFN γ in serum, plasma and other body fluids, as well as in cell-culture supernatant.

II. PRINCIPLE OF THE TEST

The Pelikine compact™ human IFN γ ELISA kit is a "sandwich-type" of enzyme immunoassay in which a monoclonal anti-human IFN γ antibody is bound onto polystyrene microtiter wells. Human IFN γ , present in a measured volume of sample or standard is captured by the antibody on the microtiter plate, and non-bound material is removed by washing. Subsequently, a biotinylated polyclonal antibody to human IFN γ is added. This antibody binds to the IFN γ -antibody complex present in the microtiter well. Excess biotinylated antibody is removed by washing, followed by addition of horseradish peroxidase (HRP) conjugated streptavidin, which binds onto the biotinylated side of the IFN γ sandwich. After removal of non-bound HRP conjugate by washing, a substrate solution is added to the wells. A coloured product is formed in proportion to the amount of IFN γ present in the sample or standard. After the reaction has been terminated by the addition of a stop solution, absorbance is measured in a microtiter plate reader. From the absorbance of samples and those of a standard curve, the concentration of IFN γ can be determined by interpolation with the standard curve.


III. STORAGE AND STABILITY

The Pelikine compact™ human IFN γ ELISA kit should be stored at -18°C to -32°C. The performance of the kit is guaranteed until the expiration date shown on the case label.

IV. CONTENTS OF THE KIT

The PeliKine compact™ human IFN γ ELISA kit contains material sufficient for 288 tests, including standard curve samples. The reagents provided are:

Quantity	Kit component		Volume	Cap colour
1 vial	coating antibody	100-fold concentrated	375 μ l	red
1 vial	blocking reagent	50-fold concentrated	2 ml	transparent
2 vials	IFN γ standard (lyophilized)	4400 pg/ml	500 μ l	black
1 vial	biotinylated IFN γ antibody	100-fold concentrated	375 μ l	yellow
1 vial	streptavidin-HRP conjugate	10,000-fold concentrated	20 μ l	brown
1 bottle	dilution buffer	5-fold concentrated	60 ml	-
3 pcs	microtiter plates + lid	-	-	-
10pcs	plate seals	-	-	-



Typical standard curve for the PeliKine compact™ human IFN γ ELISA kit

V. PRECAUTIONS FOR USE

- 1) The PeliKine compact™ human IFN γ ELISA kit is intended *for research purposes only*.
- 2) Only use the reagents and microtiter plates supplied with the kit, do not mix reagents from different kit lots.
- 3) Handle all plasma and serum samples with care to prevent transmission of blood-borne infections.
- 4) Sodium azide inactivates HRP, do not use sodium azide-containing solutions, nor add sodium azide to the supplied materials.
- 5) All reagents contain thiomersal (0.001 % w/v) and may be toxic upon ingestion, inhalation or skin contact. Avoid contact of skin, eyes or clothing with dilution, washing or substrate buffer. In case of contact, wash skin or eyes with water and consult a physician.
- 6) The IFN γ standard contains human serum which has been found to be non-reactive for Hepatitis B surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus (HIV). Nevertheless the standard should be handled as potentially hazardous and capable of transmitting diseases.
- 7) Centrifuge all vials before use (1 minute 3000 x g).
- 8) With the exception of the substrate blank wells, do not allow wells to stand uncovered or dry for extended periods between incubation steps.

	STATIC INCUBATION	SHAKEN INCUBATION
	Calculated mean absorbance at 450 nm	
substrate blank	0	0
0 pg/ml	0.025	0.073
2.0 pg/ml	0.035	0.074
5.1 pg/ml	0.040	0.095
12.8 pg/ml	0.071	0.154
32 pg/ml	0.142	0.336
80 pg/ml	0.327	0.671
200 pg/ml	0.753	1.572
500 pg/ml	1.604	> 3.000

DO NOT USE THESE DATA TO CONSTRUCT A STANDARD CURVE FOR SAMPLE VALUE CALCULATIONS

Add 100 μ l of stop solution to all wells.

After stopping the colour is stable for maximally 30 minutes.

15. PLATE READ-OUT

Read at 450 nm in an ELISA reader.

IX. RESULTS

Substrate blank

- Record the absorbance at 450 nm for the substrate blank wells and average the duplicate values.

Standard curve

- Record the absorbance at 450 nm for each well containing standard and average the duplicate values.
- Calculate the net average absorbances by subtracting the average of the substrate blank wells.
- Plot the net average absorbances (ordinate) versus the IFN γ concentration in pg/ml (abscissa) on log-linear paper and draw the best fitting curve. An example of a standard curve is given on the next page.

Samples

- Record the absorbance at 450 nm for each standard well, and average the duplicate values.
- Calculate the net average absorbances by subtracting the average of the substrate blank wells.
- Locate the net average absorbance value found for each sample on the vertical axis and follow a horizontal line intersecting the standard curve. At the point of intersection, read the IFN γ concentration (pg/ml) from the horizontal axis. Multiply the obtained IFN γ concentration with the dilution factor of the sample and record this figure.

X. INCREASED SENSITIVITY

The assay sensitivity can be increased by a small adaption of the incubation methodology. Just follow all the instructions as stated in the assay procedure (chapter VIII), but incubate at room temperature (18-25°C) on a horizontal plate shaker at **700 \pm 100 rpm**. All incubations, with exception of the enzymatic colour development, have to be completed on the shaker. This will result in an increase in assay sensitivity, with small effects on the background levels (see figure next page).

VI. ADDITIONAL BUFFERS & SOLUTIONS REQUIRED

Coating buffer: 0.1 M Carbonate/bicarbonate buffer pH 9.6

Solution A: 1.24 g Na₂CO₃.H₂O in 100 ml distilled water

Solution B: 1.68 g NaHCO₃ in 200 ml distilled water

Take 70 ml of solution A, and add solution B until the pH is 9.6 (approximately 175 ml of solution B required)

The prepared buffer can be stored up to one week at 2-8°C.

PBS stock solution [20 x]: 0.2 M Phosphate Buffered Saline (PBS)

Dissolve	32	g	Na ₂ HPO ₄ .2H ₂ O
	6	g	NaH ₂ PO ₄ .2H ₂ O
	164	g	NaCl
	in 900 ml distilled water		

(intensive stirring and some heating will speed dissolution).

Bring the temperature of the solution back to room temperature (18-25°C) and check pH; if necessary adjust pH to 6.8 - 6.9 with concentrated HCl or NaOH, and add distilled water to a volume of 1 liter (when diluted 20 times the obtained buffer will have a pH of 7.2 - 7.4).

Add 20 mg thiomersal as preservative. Do not use sodium azide (NaN₃) since this preservative reduces the quality of the enzymatic label.

The prepared buffer can be stored up to three months at 2-8°C.

Note: in the concentrated buffer salt crystals may appear when stored at 2-8°C. Before preparing the working-strength buffer, first warm the concentrated buffer BRIEFLY to 37°C to dissolve the precipitate.

Washing buffer: PBS with 0.005 % TWEEN 20

Make 1 liter of working-strength PBS by diluting the PBS stock solution (see above) 20-fold with distilled water.

Add 50 μ l TWEEN 20.

The prepared buffer can be stored up to one month at 2-8°C.

Substrate buffer: 0.11 M acetate buffer pH 5.5

Dissolve 15.0 g sodium-acetate (CH₃COONa.3H₂O) in 800 ml distilled water.

Adjust pH to 5.5 with glacial acetic acid, add distilled water to a volume of 1 liter.

Do not add any preservative (e.g. merthiolate, sodium azide) since this may affect the quality of the enzymatic colour development.

The prepared buffer can be stored up to two weeks at 2-8°C

3,5,3',5'-tetramethylbenzidine (TMB) stock solution: 6 mg/ml TMB in DMSO

Dissolve 30 mg 3,5,3',5'-tetramethylbenzidine (TMB) in 5 ml dimethylsulfoxide (DMSO).

The prepared stock solution can be stored up to 1 month at **room temperature (18-25°C)** and if **protected against light**.

Hydrogen peroxide stock solution: 3% H₂O₂ solution in distilled water.

The prepared stock solution can be stored up to one month at 2-8°C.

SUBSTRATE SOLUTION

For each plate mix the following reagents:

- 12 ml substrate buffer
- 200 μ l TMB stock solution
- 12 μ l H₂O₂ stock solution

The substrate solution should be prepared just before use and has to be at room temperature (18-25°C) for optimal reproducible results.

Stop solution: 1.8 M H₂SO₄ solution in distilled water.**VII. ADDITIONAL INFORMATION****Additional materials required**

- Pipetting devices for accurate delivery of 1-10 μ l, 50 μ l, 100 μ l and 1 ml volumes.
- Beakers, flasks, cylinders necessary for preparation of reagents.
- Device for delivery of washing buffer (wash bottle / automated plate washer).
- Microtiter plate reader.

The kit contains one brown capped vial of concentrated streptavidin-poly-HRP conjugate, which must be stored at -18°C to -32°C to maintain maximal stability. The contents of the vial will not be frozen at this temperature.

Per microtiter plate, add 3 μ l streptavidin-poly-HRP conjugate to 30 ml of working-strength kit buffer just before use. **Do not prepare in advance of assay.**

Leaving the substrate blank wells empty, add 100 μ l of streptavidin-poly-HRP conjugate to all wells.

Cover the microtiter plate(s) with adhesive seal, gently agitate the microtiter plate by tapping the edge of the frame for a few seconds to mix contents of each well and **incubate for 30 minutes at room temperature (18-25°C)**.

12. FOURTH WASH STEP

Aspirate supernatant from wells and wash the microtiter plate(s) as described in point 6.

13. FOURTH INCUBATION STEP**Enzymatic colour development**

Approximately 10 minutes before use, prepare the substrate solution as described on page 4 of this leaflet.

The substrate solution should be at room temperature (18-25°C) for optimal reproducible results.

Add 100 μ l of substrate solution to all wells, **including the substrate blank wells**.

Cover microtiter plate(s) with lid, gently agitate the microtiter plate by tapping the edge of the frame for a few seconds to mix contents of each well and **incubate for 30 minutes at room temperature (18-25°C) in the dark**.

do not cover the plate with aluminium foil.

Note: The speed of enzymatic colour development is influenced by many factors including temperature and quality of the used TMB and H₂O₂.

14. STOP ENZYMATIC REACTION

Prepare washing buffer as described on page 3 of this leaflet.

Wash the required microtiter plates five times with washing buffer in a plate washer. In case of manual washing, completely fill the wells (> 300 μ l) with washing buffer and aspirate, repeat this four times. After the final aspiration the wells should be dry.

7. FIRST INCUBATION STEP

Standards and samples

Leaving the substrate blank wells empty, transfer 100 μ l of the prepared standards, control and samples in duplicate into the appropriate wells (see recommended plate plan). Cover plate(s) with adhesive seal, gently agitate the microtiter plate by tapping the edge of the plate for a few seconds to mix contents of each well and **incubate for 1 hour at room temperature (18-25°C)**.

8. SECOND WASH STEP

Aspirate supernatant from wells and wash the microtiter plate(s) as described in point 6 above.

9. SECOND INCUBATION STEP

biotinylated antibody

The kit contains one yellow-capped vial with concentrated antibody-biotin conjugate.

Per microtiter plate, add 120 μ l biotinylated antibody to 12 ml working-strength kit buffer just before use.

Leaving the substrate blank wells empty, add 100 μ l of diluted biotinylated antibody to all wells.

Cover plate(s) with adhesive seal, gently agitate the microtiter plate by tapping the edge of the plate for a few seconds to mix contents of each well and **incubate for 1 hour at room temperature (18-25°C)**.

10. THIRD WASH STEP

Aspirate supernatant from wells and wash the microtiter plate(s) as described in point 6.

11. THIRD INCUBATION STEP

Streptavidin-poly-HRP conjugate

Sensitivity

MEAN calculated zero signal + 3 SD : 1-2 pg/ml (shake - static incubation)
2 x (MEAN calculated zero signal) : 4-6 pg/ml (shake - static incubation)

Note: the sensitivity is dependent of the type and quality of enzymatic substrate.

Expected values

IFN γ values in fresh serum and plasma samples of healthy individuals are below 10 pg/ml.

Specificity

No crossreactivity was observed with the following recombinant human proteins: IL-1 α , IL-1 β , IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-11, Macrophage Colony Stimulating Factor (M-CSF), Granulocyte Colony Stimulating Factor (G-CSF), Granulocyte/Macrophage Colony Stimulating Factor (GM-CSF), Leukaemia Inhibitory Factor (LIF), RANTES, Stem Cell Factor/Mast Cell Factor (SCF/MCF), Transforming Growth Factor β -1 (TGF β -1), Tumour Necrosis Factor α (TNF α) and Tumour Necrosis Factor β (TNF β /Lymphotoxin).

REFERENCES

- Adolf, G.R. (1985) *Oncology* **42**: 33
 Aggerwal, B.B. Gutterman, J.U. (1992) *Human cytokines*. Blackwell Sci.Pub. ISBN 0-86542-183-8
 Balkwill, F. (1989) *The Lancet* (I): 1060
 Billiau, A. (1988) *Immunology Today* **9**: 37
 Billiau, A. *et al* (1990) *Biochem.Pharmacol.* **40**: 1433
 Bruserud, O. *et al* (1993) *Eur.J.Hematol.* **51**: 73
 Celada, A. *et al* (1989) *Eur.J.Immunol.* **19**: 1103
 Doldi, K. *et al* (1985) *J.Interferon Res.* **5**: 55
 Farrar, M.A. *et al* (1993) *Ann.Rev.Immunol.* **11**: 571
 Gray, P.W. *et al* (1892) *Nature* **295**: 503
 Grossberg, S.E. *et al* (1989) *Experientia* **45**: 508
 Ijzermans, J.M. *et al* (1989) *Immunobiol.* **179**: 456
 Kwiatkowski, D.A. *et al* (1990) *The Lancet* **336**: 1201
 Locksley, R.M. *et al* (1991) *Immunoparasitology Today*: A58
 Mogensen, S.C. *et al* (1987) *Interferon* **8**: 55
 O'Garra, A. (1989) *The Lancet* (I): 1003
 Paliard, X. (1988) *et al* *J.Immunol.* **141**: 849
 Reiter, Z. (1993) *J.Interferon Res.* **13**: 247
 Samuel, C.E. (1991) *Virology* **183**: 1
 Thomson, A.W. (1991) *The cytokine handbook*. Academic Press ISBN 0-12-689660-7

VIII. ASSAY PROCEDURE

1. **BRING ALL REAGENTS TO ROOM TEMPERATURE (18-25°C)**, with the exception of the

streptavidin-HRP conjugate, positive control and standard which has to be kept at -18°C to -32°C to ensure stability. Centrifuge all vials before use (1 minute 3000 x g).

For your convenience an easy-reference manual with check list and plate plan are available on the last pages of this leaflet.

2. DILUTION BUFFER

The kit contains one bottle with 5-fold concentrated dilution buffer.

For optimal assay results, dilute samples and standards in working-strength dilution buffer.

Calculate the quantity of dilution buffer required (approximately 15 ml undiluted buffer per microtiter plate) and prepare a working-strength solution by diluting the opalescent concentrated buffer 5 times in distilled water before use. Shake gently. The working-strength dilution buffer can be stored for up to one week at 2-8°C.

3. MICROTITER PLATES

Coating antibody

Coating

The kit contains three microtiter plates for 96 tests each, including the standard curve and control samples.

Prepare coating buffer as described on page 3 of the information leaflet.

Per microtiter plate, add 120 μ l of coating antibody (red-capped vial) to 12 ml coating buffer.

Add 100 μ l to all wells, cover microtiter plate(s) with lid and **incubate overnight at room temperature (18-25°C)**.

Washing procedure

Prepare washing buffer as described on page 3 of the information leaflet.

Aspirate supernatants from wells and completely fill the wells (> 300 μ l) with working-strength PBS/TWEEN and aspirate. Repeat this four times, after the final aspiration the wells should be dry.

Blocking procedure

The kit contains one transparent-capped vial with 2 ml blocking reagent.

Prepare blocking buffer by adding 500 μ l blocking reagent to 25 ml working-strength PBS (1:20 dilution of stock PBS as described on page 3 of the information leaflet).

Add 200 μ l blocking buffer to all wells, cover microtiter plate(s) with adhesive seal, gently agitate the microtiter plate by tapping the edge of the plate for a few seconds to mix contents of each well and **incubate for 1 hour at room temperature (18-25°C)**.

4. IFN γ STANDARD

Standard curve preparation

A natural human IFN γ standard has been calibrated against the WHO reference preparation

(IFN γ 88/606; National Institute for Biological Standards and Control, Potters Bar, Hertfordshire, U.K. 1 WHO Unit = 53 pg IFN γ).

The kit contains two lyophilized black-capped vials with 4400 pg/ml natural IFN γ .

Reconstitute one lyophilized standard by adding 500 μ l of distilled water to the vial. Incubate for 10 minutes at room temperature and mix gently. After reconstitution the standard must be kept cold (2-8°C) and stored frozen after use (<-18°C, preferably <-70°C).

Label 7 tubes, one tube for each dilution: 500, 200, 80, 32, 12.8, 5.1, and 2.0 pg/ml. Pipette 585 μ l of working-strength dilution buffer into the tube labelled 500 pg/ml and 300 μ l of working-strength dilution buffer into the other tubes. Transfer 75 μ l of the IFN γ standard (4400 pg/ml) into the first tube labelled 500 pg/ml, mix well and transfer 200 μ l of this dilution into the second tube labelled 200 pg/ml. Repeat the serial dilutions six more times by adding 200 μ l of the previous tube of diluted standard to the 300 μ l of dilution buffer. The standard curve will contain 500, 200, 80, 32, 12.8, 5.1, 2.0 and 0 pg/ml (dilution buffer).

It is recommended to prepare two separate series for each assay.

Avoid repeated freeze-thawing of the standard, although experimental data have shown that up to 3 freeze-thaw cycles have no effect on the IFN γ levels of the standard. Thaw the reconstituted standard in tap water (18-25°C), do not use 37°C or 56°C water baths for this purpose. The second vial of standard can be used in case of prolonged storage of the reconstituted standard (> 2 months).

5. SAMPLES

Serum, EDTA-anti-coagulated plasmas, and culture fluids are suitable for use in the assay (**caution:** separate plasma/serum and blood cells within 4 hours after collection, non-separated samples must be kept on temperatures from 2 to 8°C). Do not use grossly haemolyzed or lipemic specimens. If samples are to be run within 24 hours, they may be stored at 2-8°C; otherwise samples should be stored frozen (<-18°C, preferably <-70°C).

Avoid freezing and thawing samples more than once. Prior to the assay, frozen samples should be thawed **as quickly as possible** in tap water (18-25°C), do not use 37°C or 56°C water baths for this purpose.

It is recommended to dilute the test samples at least 1:2 in working-strength dilution buffer. If high levels of IFN γ (> 500 pg/ml) are expected in the test samples, additional dilutions of sample i.e. 1:10 and 1:50 should also be prepared.

6. FIRST WASH STEP

